

DRESS ORTHOSES

Specifications



SYSTEM 3.0 WOMEN'S FLAT/PUMP

- **Semi-flexible Shell:** Procarbolene (2.2mm)
- **Semi-rigid Shell:** Procarbolene (2.6mm)
- **Rigid Shell:** Procarbolene (3mm)
- Bunions
- Heel spurs
- Moderate FF varus
- Dress shoes under 1.5in heel height

Semi-Flexible

Patient Weight: 175
Flex: 43
Grind: Pump
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble

Semi-Rigid

Patient Weight: 220
Flex: 32
Grind: Pump
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble

Rigid

Patient Weight: 220+
Flex: 19.5
Grind: Pump
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble



SYSTEM 3.0 WOMEN'S HIGH HEEL

- **Semi-flexible Shell:** Procarbolene (2.2mm)
- **Semi-rigid Shell:** Procarbolene (2.6mm)
- General foot fatigue
- Heel spurs
- Moderate FF valgus or varus
- Dress shoes over 1.5in heel height

Semi-Flexible

Patient Weight: 175
Flex: 43
Grind: Pump
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble

Semi-Rigid

Patient Weight: 220
Flex: 32
Grind: Pump
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble



STEEPIN' OUT

- **1/8 in Shell:** Polypropylene (3.2mm)
- Plantar fasciitis
- Moderate FF valgus or varus
- Pes cavus
- Pronation

Patient Weight:

Flex: 45
Grind: Pump
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble

DRESS ORTHOSES

Specifications



SYSTEM 3.0 MEN'S

- **Semi-flexible Shell:** Procarbolen (2.2mm)
- **Semi-rigid Shell:** Procarbolen (2.6mm)
- **Rigid Shell:** Procarbolen (3mm) • Heel spurs
- General foot fatigue
- Heel spurs
- Moderate FF valgus or varus

Semi-Flexible

Patient Weight: 175
Flex: 43
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble

Semi-Rigid

Patient Weight: 220
Flex: 32
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble

Rigid

Patient Weight: 220+
Flex: 19.5
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble



ADVANTAGE+ WOMEN'S DRESS

- **1/8 in Shell:** Polypropylene (3.2mm)
- **3/16 in Shell:** Polypropylene (4.7mm)
- General foot fatigue
- Bunions
- Metatarsalgia

Advantage+ Sport (1/8in)

Patient Weight: 175
Flex: 45
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble

Advantage+ Sport (3/16in)

Patient Weight: 220+
Flex: 34
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble



ADVANTAGE+ MEN'S DRESS

- **1/8 in Shell:** Polypropylene (3.2mm)
- **3/16 in Shell:** Polypropylene (4.7mm)
- General foot fatigue
- Heel spurs
- Moderate FF valgus or varus

Advantage+ Sport (1/8in)

Patient Weight: 175
Flex: 45
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble

Advantage+ Sport (3/16in)

Patient Weight: 220+
Flex: 34
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Intrinsic
Topcover: Blue Marble