

# LEISURE ORTHOSES

## Specifications



### PEDESTRIAN

- Shell: Polyethelene (3mm)
- Rheumatoid arthritis
- Pes cavus
- Shin splints
- General foot fatigue
- DJD
- Pronation

**Patient Weight:** 185+  
**Flex:** 47  
**Grind:** Normal  
**Heel Cup Depth:** 10mm  
**Rearfoot Post:** Extrinsic  
**Topcover:** Sand Straw



### ENHANCED FIT

- Shell: Polypropylene (3.2mm)
- Heel spurs
- Pronation
- Plantar Faciitis
- High RF varus

**Patient Weight:** 175  
**Flex:** 45  
**Grind:** Normal  
**Heel Cup Depth:** 10mm  
**Rearfoot Post:** Extrinsic  
**Topcover:** Sand Straw