

PRE-FABRICATED ORTHOSES

Specifications



TRANSITION SPORT

- Shell: thin, composite, 3/4-length Implex
- Heel Spurs
- Tendonitis
- Moderate pronation

Rearfoot Post: Neutral
Width: Normal
Topcover: Dark red Implus foam
Women Sizes: 6-12
Men Sizes: 7-14
Other: 100% Memory recovery



TRANSITION WOMEN'S DRESS

- Shell: thin, composite, 3/4-length Implex
- Bunions
- Heel Spurs
- Moderate FF valgus or varus

Rearfoot Post: Neutral extrinsic
Width: Normal
Topcover: Dark red Implus foam
Sizes: 6-12
Other: 100% Memory recovery



TRANSITION MEN'S DRESS

- Shell: thin, composite, 3/4-length Implex
- General foot fatigue
- Heel Spurs
- Moderate FF valgus or varus

Rearfoot Post: None
Width: Normal
Topcover: Dark red Implus foam
Sizes: 7-14
Other: 100% Memory recovery



CASUAL DRESS

- Polypropylene Shell
- Men Sizes: 4-13
(half sizes available excluding 12-13)
- Women Sizes: 6-13.5
(half sizes available)



PERFORMANCE RX

- Polypropylene Shell
- Men Sizes: 4-13
(half sizes available excluding 12-13)
- Women Sizes: 6-13.5
(half sizes available)

PRE-FABRICATED ORTHOSES

Specifications



FULL-LENGTH ORTHOSES

- Polypropylene Shell
- Men Sizes: 4–15
(half sizes available excluding 12-15)
- Women Sizes: 6–13.5
(half sizes available)



3/4 LENGTH ORTHOSES

- Polypropylene Shell
- Men Sizes: 3–14
(half sizes available excluding 12-14)
- Women Sizes: 5–13.5
(half sizes available)



HEAT MOLDABLE

- EVA Shell
- Men Sizes: 4–15
(half sizes available excluding 12-15)
- Women Sizes: 6–13.5
(half sizes available)



POLYPRO SHELL

- EVA Shell
- Men Sizes: 4–15
(half sizes available excluding 12-15)
- Women Sizes: 6–13.5
(half sizes available)