

SPORT ORTHOSES

Specifications



MARATHOTIC

- Shell: Polyethylene (5mm)
- Moderate FF valgus or varus
- Chondromalacia
- Postural lower back problems

Patient Weight: 180+
Flex: 39
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble



TL-ACTION

- Shell: TL-2100 Graphite (2.2mm)
- Severe FF varus
- Tarsal tunnel syndrome
- Equinus

Patient Weight: 190
Flex: 20
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble



TL-TUFF

- Shell: TL-2100 Graphite (2.7mm)
- Tarsal coalition
- Severe FF varus
- Calcaneal varus

Patient Weight: 190+
Flex: 19.5
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble



PRO-TECH

- 1/8 in Shell: Polypropylene (3.2mm)
- 3/16 in Shell: Polypropylene (4.7mm)
- Moderate FF valgus or varus
- Chondromalacia
- Postural lower back problems

1/8in
Patient Weight: 175
Flex: 45
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble

3/16in
Patient Weight: 220+
Flex: 34
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble

SPORT ORTHOSES

Specifications



SYSTEM 3.0

- **Semi-flexible Shell:** Procarbolene (2.2mm)
- **Semi-rigid Shell:** Procarbolene (2.6mm)
- **Rigid Shell:** Procarbolene (3mm) • Heel spurs
- Tendonitis
- Moderate pronation
- Shin Splints
- Plantar Fasciitis

Semi-Flexible

Patient Weight: 175
Flex: 43
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble

Semi-Rigid

Patient Weight: 220
Flex: 32
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble

Rigid

Patient Weight: 220+
Flex: 19.5
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble



ADVANTAGE+ SPORT

- **1/8 in Shell:** Polypropylene (3.2mm)
- **3/16 in Shell:** Polypropylene (4.7mm)
- Plantar fasciitis
- Heel spurs
- Shin splints
- Pronation

1/8in

Patient Weight: 175
Flex: 45
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble

3/16in

Patient Weight: 220+
Flex: 34
Grind: Normal
Heel Cup Depth: 10mm
Rearfoot Post: Extrinsic
Topcover: Burgundy Marble