


XFIT DIASYSTEM INSERTS

Specifications



XFIT DIASYSTEM INSERT

- Shell: EVA
- Shell Durometer: 35
- Neuropathy
- Toe amputations
- Metatarsalgia
- Heel pain
- Rheumatoid arthritis
- Fat atrophy

Topcover: 1/8" Plastazote
Topcover Durometer: 20
Topcover Color: Pink
Middle Layer: 1/8" PPT
Middle Layer Durometer: 20
Middle Layer Color: Blue
Patient Weight: Less than 200 lbs.



XFIT DIASYSTEM PLUS INSERT

- Shell: Cork
- Shell Durometer: 60
- Neuropathy
- Toe amputations
- Metatarsalgia
- Heel pain
- Rheumatoid arthritis
- Fat atrophy

Topcover: 1/8" Plastazote
Topcover Durometer: 20
Topcover Color: Pink
Middle Layer: 1/8" PPT
Middle Layer Durometer: 20
Middle Layer Color: Blue
Patient Weight: 200+ lbs
Other: 3 layers